

**Trichoplax *in-situ* hybridization protocol**

DATE: \_\_\_\_\_

Hybe Temp: \_\_\_\_\_

Condition	Anti-sense probe	Sense probe	Other control

- Fixation**
- Transfer animals from slides to gelatin coated dishes
  - Allow animals to settle on the bottom of the dish (at least 1 hour)
  - Fix animals by gently adding ice cold fix (4% PFA, 0.2% Glut, in high salt seawater 0.5g NaCl / 50ml) to the dish for 90 seconds.
  - Gently remove solution, add ice cold 4% PFA (in high salt seawater), place at 4 C for 1 hour.
  - Remove Fix, wash 3x Depc H2O. (ICE COLD)
  - 5 min 25% methanol-75% Depc H2O (ICE COLD)
  - 5 min 50% methanol-50% Depc H2O (ICE COLD)
  - 5 min 75% methanol-25% Depc H2O (ICE COLD)
  - Store in 100% methanol (ICE COLD)
- Pretreatment**
- Transfer animals (1 per well) to a 96 well plate(150µl / wash)
- Rehydrate**
- 5 min 75% methanol-25% PTw
  - 5 min 50% methanol-50% PTw
  - 5 min 25% methanol-75% PTw
  - Wash 5x PTw
- Proteinase**
- 5 min Proteinase K (0.01 mg/ml) at RT [5 µl per 10 ml]
  - 2x 5 min PTw + 2 mg/ml glycine at RT [100 mg/ml stock; 500 µl per 25 ml]
  - 2x 5 min 0.1% triethanolamine in water [332 µl per 25 ml]
  - Premix 0.1% triethanolamine with 3 µl/ml acetic anhydride and add to well for 5 min
  - Repeat previous step but with 6 µl/ml acetic anhydride
  - 2x 5 min PTw
- Refix**
- 1 hour 4% paraformaldehyde in PTw at 4 C [1 ml per 4 ml]
  - Wash 5x in PTw
- Prehybe**
- 2 times 10 min hybe buffer
  - Hybe buffer 3 hrs to overnight at hybe temp ⊗ [save this and previous used hybe buffer for washes]
- Hybe**
- Add probe (1 ng/µl) and hybridize overnight or the weekend ⊗
- Washes**
- Remove probe
  - 10 min 100% hybe buffer [used] at hybe temp
  - 20 min 100% hybe buffer [used] at hybe temp
  - 20 min 75% hybe - 25% 2x SSC at hybe temp
  - 20 min 50% hybe - 50% 2x SSC at hybe temp

- 20 min 25% hybe - 75% 2x SSC at hybe temp
- 3x 20 min 2x SSC at hybe temp
- 3x 10 min 0.05x SSC at hybe temp
- 5 min 75% 0.05x SSC 25% PTw at room temp
- 5 min 50% - 50% PTw
- 5 min 25% - 75% PTw
- 3x 10 min PTw
- Block** - 1 hour (or longer) Blocking Buffer at RT
- Antibody** - Incubate with anti-Dig/AP (1:5000) at overnight at 4°C ⊗
- Washes** - 10x 10 (to 30) minutes PTw or more  
3-5x 5 min PBS only
- Develop** - 2x 5 min AP Buffer WITHOUT MgCl<sub>2</sub> (USE FRESH TWEEN)
- 2x 5 min AP Buffer (USE FRESH TWEEN) – Swirl and change quickly
- Stain in AP Buffer with NBT (3.3 µl/ml) and BCIP (3.3 µl/ml)
- To STOP, wash with PBS or PTw 3-5 times
- Mount in 70% glycerol in PBS (can go through 30% glycerol first)

**Hybridization Buffer (Total: 40 ml)**

Formamide	20 ml	50%
20x SSC pH4.5	10 ml	5x
Heparin 20 mg/ml	0.1 ml	50 µg/ml
*20% Tween-20	0.5 ml	0.5%
20% SDS	2.0 ml	1.0%
SS DNA 10mg/ml	0.2 ml	100 µg/ml
Distilled Water	7.5 ml	

PTw – PBS + 0.1% Tween-20 (5 ml 20% Tween-20 per 1 L of PBS)

PBT – PBS + 0.2% TritonX-100 + 0.1% BSA

PTr – PBS + 0.2% TritonX-100

AP Buffer – 100 mM NaCl + 50 mM MgCl<sub>2</sub> + 100 mM Tris pH9.5 + 0.5% Tween-20

AP minus Mg – 100 mM NaCl + 100 mM Tris pH 9.5 + 0.5% Tween-20

**AP Stock Solution (Total: 10 ml)**

Distilled Water	7.25 ml
1M NaCl	1.00 ml
1M Tris pH9.5	1.00 ml
20% Tween-20	0.25 ml

1M MgCl<sub>2</sub> (or water) 0.50 ml

MAKE UP FRESH 20% Tween-20!